

MOCS3 Antibody

Catalog No: #40000


Package Size: #40000 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MOCS3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Species Reactivity	Human;Mouse;Rat
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal MOCS3 (1) antibody was raised against a recombinant human MOCS3 protein 15-175aa (BC015939).
Conjugates	Unconjugated
Target Name	MOCS3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% sodium azide.
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.