

mOrange Monoclonal Antibody

Catalog No: #40480



Package Size: #40480-1 50ul #40480-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	mOrange Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Specificity	The antibody detects mOrange & mOrange Tag protein.
Target Name	mOrange
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:5000

Images



Western blot analysis of recombinant mOrange protein,
diluted at 1:10000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.