

Annexin II Polyclonal Antibody

Catalog No: #40592



Package Size: #40592-1 50ul #40592-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin II Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Annexin II Polyclonal Antibody detects endogenous levels of Annexin II protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Annexin II.
Conjugates	Unconjugated
Target Name	Annexin II
Other Names	ANXA2; ANX2; ANX2L4; CAL1H; LPC2D; Annexin A2; Annexin II; Annexin-2; Calpactin I heavy chain; Calpactin-1 heavy chain; Chromobindin-8; Lipocortin II; Placental anticoagulant protein IV; PAP-IV; Protein I; p36
Accession No.	Swiss-Prot: P07355NCBI Gene ID: 302
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Annexin II Polyclonal Antibody



Western Blot analysis of Jurkat NIH-3T3 cells using Annexin II Polyclonal Antibody

Published Papers

el at., LncRNA-MEG3 attenuates hyperglycemia-induced damage by enhancing mitochondrial translocation of HSP90A in the primary hippocampal neurons. In Exp Cell Res on 2022 Oct 15 by Wanqing Zhou, Zhihua Wang,et al..PMID:35998683, , (2022)
[PMID:35998683](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.