

NFkB-p100 Polyclonal Antibody

Catalog No: #41224



Package Size: #41224-1 50ul #41224-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NFkB-p100 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IP;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	NFkB-p100 Polyclonal Antibody detects endogenous levels of NFkB-p100 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human NFkB-p100 around the non-phosphorylation site of S869.
Conjugates	Unconjugated
Target Name	NFkB-p100
Other Names	NFKB2; LYT10; Nuclear factor NF-kappa-B p100 subunit; DNA-binding factor KBF2; H2TF1; Lymphocyte translocation chromosome 10 protein; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2; Oncogene Lym-10; Lym10
Accession No.	Swiss-Prot: Q00653NCBI Gene ID: 4791
Calculated MW	97kDa
SDS-PAGE MW	110-120kDa
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-1:2000; IHC 1:100-1:300; IF 1:50-300; IP 2-5ug/mg lysate; ELISA 1:20000;

Images



Western Blot analysis of Jurkat cells using NFkB-p100 Polyclonal Antibody



Western Blot analysis of Jurkat cells using NFκB-p100 Polyclonal Antibody

Published Papers

el at., Single-cell analysis reveals the immune heterogeneity and interactions in lungs undergoing hepatic ischemia-reperfusionInInt ImmunopharmacolOn2023 NovbyChen Zhang1,Mingwei Sheng2, et al..PMID:37844464, , (2023)
[PMID:37844464](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.