

BRMS-1 Polyclonal Antibody

Catalog No: #41571

Package Size: #41571-1 50ul #41571-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BRMS-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	BRMS-1 Polyclonal Antibody detects endogenous levels of BRMS-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human BRMS-1.
Conjugates	Unconjugated
Target Name	BRMS-1
Other Names	BRMS1; Breast cancer metastasis-suppressor 1
Accession No.	Swiss-Prot: Q9HCU9NCBI Gene ID: 25855
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 K562 BT474 cells using BRMS-1 Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.