

FGA Polyclonal Antibody

Catalog No: #42168

Package Size: #42168 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FGA Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total FGA polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human FGA proteinB£B"36-218aaB£B©
Target Name	FGA
Other Names	FGA protein
Accession No.	Swiss-Prot#: Q8WW76
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20 - 1:200

Images



Immunohistochemical analysis of paraffin-embedded human placenta using #42168 at dilution of 1:100.

References

[1]The status, quality, and expansion of the NIH full-length cDNA project: the Mammalian Gene Collection (MGC). The MGC Project Team Gerhard D.S., Wagner L., Feingold E.A., Shenmen C.M., Grouse L.H., Schuler G., Klein S.L., Old S., Rasooly R., Good P., Guyer M., Peck A.M., Derge J.G., Lipman D., Collins F.S., Jang W., Sherry S., Feolo M. expand/collapse author list Malek J. Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.