

GFP-Tag Monoclonal Antibody(1G6)

Catalog No: #44011



Package Size: #44011-1 50ul #44011-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GFP-Tag Monoclonal Antibody(1G6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1G6
Isotype	IgG1
Applications	WB
Specificity	The antibody detects GFP and GFP tag fusion proteins.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	GFP-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



2ug GFP fusion protein+ Primary antibody dilution at
1) 1:5000
2) 1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.