

TNF α Monoclonal antibody(Q34)

Catalog No: #44073

Package Size: #44073-1 50ul #44073-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF α Monoclonal antibody(Q34)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Q34
Isotype	IgG1
Applications	WB IHC IF
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous TNF α protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	TNF- α
Other Names	TNF- α
Accession No.	P01375
Calculated MW	17,25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-1:3000

Immunohistochemistry: 1:50-1:200

IF 1:50-200

Images



Western blot analysis of Recombinant Human TNF α Protein with TNF α Mouse mAb at 1:2,000.


The picture was kindly provided by our customer

mouse BV2



Huazhong University of Science and Technology (HUST)

Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using TNF α Mouse mAb diluted at 1:50



Published Papers

el at., Screening of anti-ζ Ahronic nonbacterial prostatitis activity of different extractions of the aerial part of Glycyrrhiza uralensis, and network pharmacology research. In Biomed Rep on 2021 Dec by Haifan Liu, Jie Cui,et al..PMID:34667596, , (2021)

[PMID:34667596](#)

el at., Screening of anti-ζ Ahronic nonbacterial prostatitis activity of different extractions of the aerial part of Glycyrrhiza uralensis, and network pharmacology research. In Biomed Rep on 2021 Dec by Haifan Liu, Jie Cui,et al..PMID:34667596, , (2021)

[PMID:34667596](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.