

IL-8 Mouse Monoclonal Antibody(14E1)

Catalog No: #44133

Package Size: #44133-1 50ul #44133-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL-8 Mouse Monoclonal Antibody(14E1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1.40E+02
Isotype	IgG1
Applications	IHC
Species Reactivity	Human;Mouse;Rat
Specificity	IL-8 Mouse Monoclonal Antibody(14E1) detects endogenous levels of IL8
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	IL-8
Other Names	IL8
Calculated MW	11kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:50-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using IL-8 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using IL-8 Mouse mAb diluted at 1:200.

Published Papers

el at., The Crosstalk Between Cancer Cells and Neutrophils Enhances Hepatocellular Carcinoma Metastasis via Neutrophil Extracellular Traps-Associated Cathepsin G Component:

A Potential Therapeutic Target. In J Hepatocell Carcinoma on 2021 May 20 by Xiangqian Guan, Yuyan Lu,et al..PMID: 34046369, , (2021)

[PMID:34046369](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.