

# Bax Mouse Monoclonal Antibody(1C1)

Catalog No: #44150

Package Size: #44150-1 50ul #44150-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                         |
|-----------------------|-------------------------------------------------------------------------|
| Product Name          | Bax Mouse Monoclonal Antibody(1C1)                                      |
| Host Species          | Mouse                                                                   |
| Clonality             | Monoclonal                                                              |
| Clone No.             | 1C1                                                                     |
| Isotype               | IgG1                                                                    |
| Applications          | IHC                                                                     |
| Species Reactivity    | Human;Mouse;Rat                                                         |
| Specificity           | Bax Mouse Monoclonal Antibody(1C1) detects endogenous levels of BAX     |
| Immunogen Type        | Peptide                                                                 |
| Immunogen Description | Synthetic Peptide                                                       |
| Conjugates            | Unconjugated                                                            |
| Target Name           | Bax                                                                     |
| Other Names           | BAX                                                                     |
| Calculated MW         | 20kd                                                                    |
| Concentration         | 1mg/ml                                                                  |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage               | Store at -20°C/1 year                                                   |


## Application Details

Immunohistochemistry: 1:100-1:200


## Images



Immunohistochemical analysis of paraffin-embedded Rat Testis Tissue using Bax Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using Bax Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Bax Mouse mAb diluted at 1:200.

## Published Papers

JuBo Wang;DanFeng Dong;MaoDe Wang;Ke Gao et al., IDH1 Overexpression Induced Chemotherapy Resistance and IDH1 Mutation Enhanced Chemotherapy Sensitivity in Glioma Cells in Vitro and in Vivo, , (2014)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.