

## TMC8 Antibody

Catalog No: #44426



Package Size: #44426-1 50ul #44426-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TMC8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	TMC8 Antibody detects endogenous levels of total TMC8
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TMC8
Conjugates	Unconjugated
Target Name	TMC8
Other Names	TMC8, EVER2, EVIN2, Transmembrane channel-like 8, EV2
Accession No.	Swiss-Prot#: Q8IU68NCBI Gene ID:147138
Calculated MW	81kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

## Published Papers

et al., Transmembrane channel-like protein 8 as a potential biomarker for poor prognosis of hepatocellular carcinoma. In Mol Clin Oncol on 2017 Aug by Panpan Lu, Qiang Ding, et al.. PMID: 28781795, (2017)  
[PMID:28781795](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.