

## IL11 Antibody

Catalog No: #44814



Package Size: #44814-1 50ul #44814-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | IL11 Antibody   |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.   |
| Applications          | WB  |
| Species Reactivity    | Human;Mouse;Rat   |
| Specificity           | IL11 antibody detects endogenous levels of total IL11   |
| Immunogen Type        | Peptide   |
| Immunogen Description | A synthesized peptide derived from human IL11   |
| Conjugates            | Unconjugated  |
| Target Name           | IL11  |
| Other Names           | Adipogenesis inhibitory factor antibody<br>AGIF antibody<br>IL 11 antibody<br>IL-11 antibody<br>IL11 antibody<br>IL11_HUMAN antibody<br>Interleukin 11 antibody<br>Interleukin-11 antibody<br>Oprelvekin antibody |
| Accession No.         | Swiss-Prot#: P20809NCBI Gene ID:3589  |
| Calculated MW         | 27 kDa?   |
| Concentration         | 1.0mg/mL  |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.  |
| Storage               | Store at -20°C  |

## Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.