

Perilipin A Antibody

Catalog No: #45080

Package Size: #45080-1 50ul #45080-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Perilipin A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	Perilipin A Antibody detects endogenous levels of total Perilipin A
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Perilipin A
Conjugates	Unconjugated
Target Name	Perilipin A
Other Names	PLIN1, FPLD4, PERI, Perilipin-1, PLIN, Perilipin, Perilipin 1
Accession No.	Swiss-Prot#: O60240NCBI Gene ID:5346
Calculated MW	62 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Published Papers

el at., FTY720 Improves the Survival of Autologous Fat Grafting by Modulating Macrophages Toward M2 Polarization Via STAT3 Pathway In Cell Transplant on Jan-Dec 2021 by Yi Yi, Weijie Hu,et al..PMID:34662222, , (2021)

[PMID:34662222](https://pubmed.ncbi.nlm.nih.gov/34662222/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.