

## GCNT4 Antibody

Catalog No: #45802



Package Size: #45802-1 50ul #45802-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GCNT4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	GCNT4 Antibody detects endogenous levels of total GCNT4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GCNT4
Conjugates	Unconjugated
Target Name	GCNT4
Other Names	GCNT4, Core2-GlcNAc-transferase 3, Core 2-branching enzyme 3, C2GNT3
Accession No.	Swiss-Prot#: Q9P109NCBI Gene ID:51301
Calculated MW	53kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis GCNT4 using HuvEc whole cell lysates

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.