

GPSM3 Antibody

Catalog No: #45995

Package Size: #45995-1 50ul #45995-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPSM3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	GPSM3 Antibody detects endogenous levels of total GPSM3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPSM3
Conjugates	Unconjugated
Target Name	GPSM3
Other Names	AGS4, C6orf9, G18
Accession No.	Swiss-Prot#: Q9Y4H4NCBI Gene ID:63940
Calculated MW	18 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis GPSM3 using COLO205 whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.