

CASR Antibody

Catalog No: #46408

Package Size: #46408-1 50ul #46408-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CASR Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CASR protein.
Immunogen Type	peptide
Immunogen Description	Synthetic protein corresponding to residues near the C terminal of human CASR
Target Name	CASR
Other Names	CAR; FHH; FIH; HHC; EIG8; HHC1; NSHPT; PCAR1; GPRC2A; HYPOC1
Accession No.	Swiss-Prot:P41180NCBI Gene ID:846NCBI Protein:BC112236
Calculated MW	121 kDa
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000


Immunohistochemistry: 1: 50-300

Images



Gel: 6%SDS-PAGE

lysate: 40 B₁ g, Lane: Hela cell lysate,Primary antibody: 46408B₁CASR Antibody) at dilution 1/450Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,
Exposure time: 2 minutes



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46408(CASR Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: x200)

Background

The protein encoded by this gene is a G protein-coupled receptor that is expressed in the parathyroid hormone (PTH)-producing chief cells of the parathyroid gland, and the cells lining the kidney tubule. It senses small changes in circulating calcium concentration and couples this information to intracellular signaling pathways that modify PTH secretion or renal cation handling, thus this protein plays an essential role in maintaining mineral ion homeostasis. Mutations in this gene cause familial hypocalciuric hypercalcemia, familial, isolated hypoparathyroidism, and neonatal severe primary hyperparathyroidism.

Published Papers

el at., Acetylcholine protects mesenteric arteries against hypoxia/reoxygenation injury via inhibiting calcium-sensing receptor.In J Pharmacol Sci on 2015 Apr by Ming Zhao , Xi He et al..PMID:25922231 , , (2015)

[PMID:25922231](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.