

CCDC137 Antibody

Catalog No: #46419



Package Size: #46419-1 50ul #46419-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CCDC137 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CCDC137 protein.
Immunogen Type	peptide
Immunogen Description	Full length fusion protein of human CCDC137
Target Name	CCDC137
Other Names	RaRF
Accession No.	Swiss-Prot:Q6PK04NCBI Gene ID:339230NCBI Protein:BC009369
Calculated MW	33 kDa
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Immunohistochemistry: 1: 30-150

Images




Gel: 8%SDS-PAGE

lysate: 40 B|Γ g, Lane: 293T cell lysate,

Primary antibody: 46419B£B`CCDC137 Antibody) at dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 46419(CCDC137 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: x200)

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.