

## GPANK1 Antibody

Catalog No: #46566

Package Size: #46566-1 50ul #46566-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GPANK1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GPANK1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic protein corresponding to residues near the C terminal of human GPANK1
Target Name	GPANK1
Other Names	G5; BAT4; D6S54E; ANKRD59; GPATCH10
Accession No.	Swiss-Prot:O95872NCBI Gene ID:7918NCBI Protein:BC008783
Concentration	0.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1: 25-100

## Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46566(GPANK1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)

## Background

This gene is located in a cluster of HLA-B-associated transcripts, which is included in the human major histocompatibility complex III region. This gene encodes a protein which is thought to play a role in immunity. Multiple alternatively spliced variants, encoding the same protein, have been identified.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.