

LDH-B Antibody

Catalog No: #48009



Package Size: #48009-1 50ul #48009-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LDH-B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB, IF-Cell, FC, IHC-P
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Synthetic peptide within Human LDH-B
Conjugates	Unconjugated
Other Names	Epididymis secretory protein Li 281 antibody HEL S 281 antibody L lactate dehydrogenase B chain antibody L-lactate dehydrogenase B chain antibody Lactate dehydrogenase H chain antibody LDH B antibody LDH H antibody LDH heart subunit antibody LDH-B antibody LDH-H antibody LDHB antibody LDHB_HUMAN antibody LDHBD antibody LDHH antibody Renal carcinoma antigen NY REN 46 antibody Renal carcinoma antigen NY-REN-46 antibody TRG-5 antibody TRG5 antibody
Accession No.	Swiss-Prot#:P07195
Calculated MW	37 kDa
Formulation	1*PBS (pH7.4), 0.2% BSA, 25% Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Application Details

WB 1:500-1:1000 IF-Cell 1:200 FC 1:500-1:1000 IHC-P 1:200-1:500

Images




Western blot analysis of LDH-B on different lysates with Rabbit anti-LDH-B antibody at 1/500 dilution.
Lane 1: Hela cell lysate
Lane 2: Jurkat cell lysate
Lane 3: Daudi cell lysate
Lane 4: HL-60 cell lysate


Lysates/proteins at 10 µg/Lane.

Exposure time: 2 minutes;


12% SDS-PAGE gel.




Immunocytochemistry analysis of A549 cells labeling LDHB with Rabbit anti-LDH-B antibody at 1/200 dilution.




Immunohistochemical analysis of paraffin-embedded rat kidney tissue with Rabbit anti-LDH-B antibody at 1/200 dilution.



Immunohistochemical analysis of paraffin-embedded human kidney tissue with Rabbit anti-LDH-B antibody at 1/200 dilution.



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue with Rabbit anti-LDH-B antibody at 1/200 dilution.



Flow cytometric analysis of A549 cells labeling LDH-B.

Background

Lactate dehydrogenase (LDH) is an enzyme present in a wide variety of organisms, including plants and animals. It catalyses the interconversion of pyruvate and lactate with concomitant interconversion of NADH and NAD⁺. In medicine, LDH is often used as a marker of tissue breakdown as LDH is abundant in red blood cells and can function as a marker for hemolysis. In mammals, three types of LDH subunits (35 kDa) are encoded by the genes Ldh-A, Ldh-B, and Ldh-C. Lactate dehydrogenase B (LDH-B, heart subunit, LDH-H) is involved in the conversion of L-lactate and NAD to pyruvate and NADH and it is predominantly localized in the heart tissue. Similar to other LDH subunit, LDH-B is considered to be an important marker for germ cell tumor.

References

1. Takeno T., Li S.S.-L.; "Structure of the human lactate dehydrogenase B gene."; Biochem. J. 257:921-924(1989).
2. The MGC Project Team; "The status, quality, and expansion of the NIH full-length cDNA project: the Mammalian Gene Collection (MGC)."; Genome Res. 14:2121-2127(2004).
3. Sakai I., Sharief F.S., Pan Y.-C.E., Li S.S.-L.; "The cDNA and protein sequences of human lactate dehydrogenase B."; Biochem. J. 248:933-936(1987).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.