

MTAP Rabbit mAb

Catalog No: #52812

Package Size: #52812-1 50ul #52812-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MTAP Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S07-5A0
Isotype	IgG
Purification	Affinity-chromatography
Applications	WB;IHC;IF/ICC;IP;FC
Species Reactivity	Human;Mouse;Rat
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MTAP
Conjugates	Unconjugated
Modification	Unmodification
Other Names	BDMF; MSAP; DMSFH; LGMBF; DMSMFH; c86fus; HEL-249
Accession No.	Swiss-Prot:Q13126GeneID:4507
Calculated MW	31 kDa
SDS-PAGE MW	31 kDa
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


WB 1:1000-1:2000; IHC 1:100-1:200; IF/ICC 1:50-1:200; IP 1:20-1:50; FC 1:20-1:100

Images



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (WB at 1:1K dilution)



Background

This gene encodes an enzyme that plays a major role in polyamine metabolism and is important for the salvage of both adenine and methionine. The encoded enzyme is deficient in many cancers because this gene and the tumor suppressor p16 gene are co-deleted. Multiple alternatively spliced transcript variants have been described for this gene, but their full-length natures remain unknown. [provided by RefSeq, Jul 2008]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.