

GRSF1 Rabbit Polyclonal Antibody

Catalog No: #54507



Package Size: #54507-1 50ul #54507-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GRSF1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant protein of human GRSF1.
Conjugates	Unconjugated
Other Names	GRSF1
Accession No.	Swiss Prot:Q12849GenelD:2926
Calculated MW	53kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using GRSF1 at 1:1000 dilution.



Immunofluorescence analysis of C6 cells using GRSF1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of L929 cells using GRSF1 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

The protein encoded by this gene is a cellular protein that binds RNAs containing the G-rich element. The protein is localized in the cytoplasm, and has been shown to stimulate translation of viral mRNAs in vitro. Multiple transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.