

CEACAM8 Rabbit Polyclonal Antibody

Catalog No: #55428



Package Size: #55428-1 50ul #55428-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CEACAM8 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human CEACAM8 (NP_001807.2).
Conjugates	Unconjugated
Other Names	CEACAM8;CD66b;CD67;CGM6;NCA-95
Accession No.	Swiss Prot:P31997GenelD:1088
Calculated MW	38kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Immunofluorescence analysis of Jurkat cells using CEACAM8 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

—

Published Papers

Xin Li et al., Alpha-1 Antitrypsin Mediates Spontaneous Resolution of Acute Gouty Inflammation Via Inhibiting Caspase 3/GSDME-dependent

Note: This product is for in vitro research use only and is not intended for use in humans or animals.