

RSAD2 Rabbit Polyclonal Antibody

Catalog No: #55466

Package Size: #55466-1 50ul #55466-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RSAD2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human RSAD2 (NP_542388.2).
Conjugates	Unconjugated
Other Names	RSAD2;2510004L01Rik;cig33;cig5;vig1
Accession No.	Swiss Prot:Q8WVG1GeneID:91543
Calculated MW	42kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using RSAD2 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.