

Biotin Rabbit pAb

Catalog No: #55648



Package Size: #55648-1 50ul #55648-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Biotin Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Biotinylated protein.
Conjugates	Unconjugated
Calculated MW	25kDa/55kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of Biotinylated Goat IgG, using Biotin Antibody.

Published Papers

Wang Shengqi;Li Jing;Xu Shang;Wang Neng;Pan Bo;Yang Bowen;Zheng Yifeng;Zhang Juping;Peng Fu;Peng Cheng;Wang Zhiyu et al., Baohuoside I chemosensitises breast cancer to paclitaxel by suppressing extracellular vesicle/CXCL1 signal released from apoptotic cells, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.