

AMBP Rabbit mAb

Catalog No: #56224



Package Size: #56224-1 50ul #56224-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	AMBP Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	AMBP Antibody detects endogenous levels of total AMBP
Immunogen Description	A synthesized peptide derived from human AMBP
Conjugates	Unconjugated
Other Names	A1M; alpha 1 microglobulin bikunin; Alpha 1 microglycoprotein; AMBP; Bikunin; EDC1; HI30;
Accession No.	Uniprot:P02760
Calculated MW	39kDa
SDS-PAGE MW	39kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

IHC:1:50~1:200

Images



All lanes: AMBP Rabbit mAb at 1/1k dilution

Lane 1 : Human plasma lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 39 kDa

Observed band size: 39 kDa

Exposure time: 6 seconds




All lanes: AMBP Rabbit mAb at 1/1k dilution

Lane 1 : Wild-type HeLa cell lysate

Lane 2 : AMBP knockdown HeLa cell lysate

Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human liver tissue stained for AMBP using 56224 at 1/100 dilution in immunohistochemical analysis.

Product Description

Inter-alpha-trypsin inhibitor inhibits trypsin, plasmin, and lysosomal granulocytic elastase. Inhibits calcium oxalate crystallization. Trypstatin is a trypsin inhibitor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.