

CD82 Rabbit mAb

Catalog No: #56256



Package Size: #56256-1 50ul #56256-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD82 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human;Rat
Specificity	CD82 Antibody detects endogenous levels of total CD82
Immunogen Description	A synthesized peptide derived from human CD82
Conjugates	Unconjugated
Other Names	C33; CD82; IA4; Inducible membrane protein; KAI1; Kangai1; SAR2; ST6; Tetraspanin 27; Tspan27;
Accession No.	Uniprot:P27701
Calculated MW	30kDa
SDS-PAGE MW	30-60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000

IHC:1:50~1:200

FC:1:50

Images



All lanes: CD82 Rabbit mAb at 1/1k dilution

Lane 1 : JK whole cell lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 30 kDa

Observed band size: 30-60 kDa

Exposure time: 7 seconds



All lanes: CD82 Rabbit mAb at 1/1k dilution

Lane 1 : Wild-type HT-1080 cell lysate

Lane 2 : CD82 knockdown HT-1080 cell lysate

Lysates/proteins at 20 µg per lane.

Product Description

Associates with CD4 or CD8 and delivers costimulatory signals for the TCR/CD3 pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.