

RhoGAP Rabbit mAb

Catalog No: #56804

Package Size: #56804-1 50ul #56804-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RhoGAP Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human;Mouse;Rat
Specificity	RhoGAP Antibody detects endogenous levels of total RhoGAP
Immunogen Description	A synthesized peptide derived from human RhoGAP
Conjugates	Unconjugated
Other Names	ARHGAP5; p190 B; p190 RhoGAP; Rho GTPase activating protein 5; RhoGAP; RhoGAP5;
Accession No.	Uniprot:Q13017
Calculated MW	172kDa
SDS-PAGE MW	172kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000

IHC:1:50~1:200

FC:1:50

Images



All lanes: RhoGAP Rabbit mAb at 1/1k dilution

Lane 1 : 3T3 whole cell lysates Lane 2 : PC12 whole cell lysates Lane 3 : HepG2 whole cell lysates Lane 4 : Hela whole cell lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 172 kDa

Observed band size: 172 kDa

Exposure time: 9 seconds




All lanes: RhoGAP Rabbit mAb at 1/1k dilution

Lane 1 : Wild-type HeLa cell lysate

Lane 2 : RhoGAP knockdown HeLa cell lysate

Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human testis tissue stained for RhoGAP using 56804 at 1/100 dilution in immunohistochemical analysis.

Product Description

GTPase-activating protein for Rho family members (PubMed:8537347).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.