

Cytokeratin 6A Rabbit mAb

Catalog No: #59006



Package Size: #59006-1 50ul #59006-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Cytokeratin 6A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC;ICC/IF
Species Reactivity	Human
Immunogen Description	A synthesized peptide derived from human Cytokeratin 6A
Conjugates	Unconjugated
Other Names	CK 6A; CK 6B; CK 6C; CK 6D; CK 6E; CK-6C; CK-6E; Cytokeratin 6a; Cytokeratin 6B; Cytokeratin 6C; Cytokeratin 6D; Cytokeratin 6E; Cytokeratin-6C; Cytokeratin-6E; K6a keratin; K6b keratin; K6C; K6c keratin; K6d keratin; K6e keratin; Keratin K6h;
Accession No.	P02538
Calculated MW	60 kDa
SDS-PAGE MW	60 kDa
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.


Application Details

WB 1:1000-1:2000IHC 1:100-1:200ICC/IF 1:50-1:200

Images



Western blot analysis of Cytokeratin 6 expression in A431 cell lysate.



Immunohistochemical analysis of paraffin-embedded human lung cancer, using Cytokeratin 6 Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.