

Sortilin Rabbit mAb

Catalog No: #59471



Package Size: #59471-1 50ul #59471-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Sortilin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Human Mouse Rat
Specificity	Sortilin Antibody detects endogenous levels of total Sortilin
Conjugates	Unconjugated
Other Names	100 kDa NT receptor; Glycoprotein 95; Gp95; LDLCQ6; Neurotensin receptor 3; NT3; NTR3; SORT 1; Sortilin 1;
Accession No.	Q99523
Calculated MW	92 kDa
SDS-PAGE MW	110 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200


Images



All lanes: Sortilin Rabbit mAb at 1/1k dilution
 Lane 1 : SW480 whole cell lysates
 Lane 2 : Hela whole cell lysates
 Lane 3 : JK whole cell lysates
 Lane 4 : C6 whole cell lysates
 Lane 5 : Mouse kidney lysates
 Lysates/proteins at 20 µg per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 92 kDa
 Observed band size: 110 kDa
 Exposure time: 7 seconds



All lanes: Sortilin Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HT-1080 cell lysate
 Lane 2 : Sortilin Rabbit mAb knockdown HT-1080 cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human brain tissue stained for Sortilin using 59471 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.