

TXNRD2 Rabbit mAb

Catalog No: #59996



Package Size: #59996-1 50ul #59996-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	TXNRD2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Human Mouse Rat
Specificity	TXNRD2 Antibody detects endogenous levels of total TXNRD2
Conjugates	Unconjugated
Other Names	SelZ; TR 3; TR; TR beta; TR3; TRXR2; TXNRD 2; Txnrd2;
Accession No.	Q9NNW7
Calculated MW	57 kDa
SDS-PAGE MW	57 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200


Images



All lanes: TXNRD2 Rabbit mAb at 1/1k dilution
Lane 1 : Hela whole cell lysates
Lane 2 : HepG2 whole cell lysates
Lane 3 : MCF-7 whole cell lysates
Lysates/proteins at 20 µg per lane.
SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 57 kDa
Observed band size: 57 kDa
Exposure time: 8 seconds



All lanes:TXNRD2 Rabbit mAb at 1/1k dilution
Lane 1 : Wild-type HT-1080 cell lysate
Lane 2 :TXNRD2 Rabbit mAb knockdown HT-1080 cell lysate
Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human kidney cancer tissue stained for TXNRD2 using 59996 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.