

TRAF6BP Rabbit mAb

Catalog No: #59997



Package Size: #59997-1 50ul #59997-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	TRAF6BP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Human Mouse Rat
Specificity	TRAF6BP Antibody detects endogenous levels of total TRAF6BP
Conjugates	Unconjugated
Other Names	CALCOCO3; D6Ert404e; D6Ert772e; PRO0105; T6BP; TAX1BP1; tax1bp1b; TXBP151;
Accession No.	Q86VP1
Calculated MW	91 kDa
SDS-PAGE MW	91 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200


Images



All lanes: TRAF6BP Rabbit mAb at 1/1k dilution
Lane 1 : HepG2 whole cell lysates
Lane 2 : Mouse testis lysates
Lane 3 : Rat testis lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 91 kDa
Observed band size: 91 kDa
Exposure time: 7 seconds



All lanes: TRAF6BP Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HAP1 cell lysate
 Lane 2 : TRAF6BP Rabbit mAb knockdown HAP1 cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human lung cancer tissue stained for TRAF6BP using 59997 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.