

PRKRA Rabbit mAb

Catalog No: #59998



Package Size: #59998-1 50ul #59998-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PRKRA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC;ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	PRKRA Antibody detects endogenous levels of total PRKRA
Conjugates	Unconjugated
Other Names	RAX; PACT; DYT16; HSD14;
Accession No.	O75569
Calculated MW	34 kDa
SDS-PAGE MW	34 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500-1:2000


IHC: 1:50-1:200

ICC/IF: 1:50-1:200


Images




All lanes: PRKRA Rabbit mAb at 1/1k dilution
Lane 1 : HepG2 whole cell lysates
Lane 2 : Mouse testis lysates
Lane 3 : Rat testis lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 34 kDa
Observed band size: 34 kDa
Exposure time: 5 seconds



All lanes:PRKRA Rabbit mAb at 1/1k dilutionLane 1 :
Wild-type HeLa cell lysateLane 2 :PRKRA Rabbit mAb
knockdown HeLa cell lysateLysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human colon cancer tissue
stained for PRKRA using 59998 at 1/100 dilution in
immunohistochemical analysis.



Immunocytochemistry/ Immunofluorescence PRKRA antibody
(59998) ICC/IF staining of PRKRA in
293T cells. Cells were fixed with 4% Paraformaldehyde
permeabilized with 0.1% Triton X-100.

Samples were
incubated with 59998 at a working dilution of 1/100. The
secondary antibody was Alexa FluorB 488 goat anti rabbit,
used at a dilution of 1/500.

Nuclei

were counterstained with DAPI.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.