

Poliovirus Receptor Rabbit mAb

Catalog No: #60005



Package Size: #60005-1 50ul #60005-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	Poliovirus Receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;ICC/IF
Species Reactivity	Human
Specificity	Poliovirus Receptor Antibody detects endogenous levels of total Poliovirus Receptor
Conjugates	Unconjugated
Other Names	CD155; HVED; mE4; Nec15; PVR; PVS; Taa1; TAGE4;
Accession No.	P15151
Calculated MW	45 kDa
SDS-PAGE MW	70 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

ICC/IF: 1:50-1:200


Images



All lanes: Poliovirus Receptor Rabbit mAb at 1/1k dilution
Lane 1 : U87-MG whole cell lysates Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 45 kDa Observed band size: 70 kDa
Exposure time: 3 seconds



All lanes: Poliovirus Receptor Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HeLa cell lysate
 Lane 2 : Poliovirus Receptor Rabbit mAb knockdown HeLa cell lysate
 Lysates/proteins at 20 µg per lane.



Immunocytochemistry/ Immunofluorescence Poliovirus Receptor antibody (60005)

ICC/IF staining of Poliovirus Receptor in HUVEC cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100.

Samples were incubated with 60005 at a working dilution of 1/100. The secondary antibody was Alexa FluorB 488 goat anti rabbit, used at a dilution of 1/500.

Nuclei were counterstained with DAPI.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.