

TEAD1 Rabbit mAb

Catalog No: #60016



Package Size: #60016-1 50ul #60016-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TEAD1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC;ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	TEAD1 Antibody detects endogenous levels of total TEAD1
Conjugates	Unconjugated
Other Names	NTEF1; Protein GT IIC; REF1; TCF13; TEA domain family member 1; TEAD1; TEF1; Transcription factor 13; Transcriptional enhancer factor TEF1;
Accession No.	P28347
Calculated MW	48 kDa
SDS-PAGE MW	52 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500-1:2000


IHC: 1:50-1:200

ICC/IF: 1:50-1:200


Images




All lanes: TEAD1 Rabbit mAb at 1/1k dilution
 Lane 1 : HeLa whole cell lysates
 Lane 2 : 293 whole cell lysates
 Lane 3 : Mouse heart lysates
 Lane 4 : Rat heart lysates
 Lysates/proteins at 20 µg per lane.
 Secondary: Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 48 kDa
 Observed band size: 52 kDa
 Exposure time: 5 seconds



All lanes: TEAD1 Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HeLa cell lysate
 Lane 2 : TEAD1 Rabbit mAb knockdown HeLa cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human colon cancer tissue stained for TEAD1 using 60016 at 1/100 dilution in immunohistochemical analysis.



Immunocytochemistry/ Immunofluorescence TEAD1 antibody (60016) ICC/IF staining of TEAD1 in HeLa cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100.

Samples were incubated with 60016 at a working dilution of 1/100. The secondary antibody was Alexa FluorB 488 goat anti rabbit, used at a dilution of 1/500.

Nuclei were counterstained with DAPI.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.