

PLAP Rabbit mAb

Catalog No: #60031



Package Size: #60031-1 50ul #60031-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PLAP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC;ICC/IF
Species Reactivity	Human
Specificity	PLAP Antibody detects endogenous levels of total PLAP
Conjugates	Unconjugated
Other Names	Glycerophosphatase; PLA2P; Plaa; PLAP;
Accession No.	Q9Y263
Calculated MW	87 kDa
SDS-PAGE MW	75 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500-1:2000


IHC: 1:50-1:200

ICC/IF: 1:50-1:200


Images




All lanes: PLAP Rabbit mAb at 1/1k dilution
 Lane 1 : MCF-7 whole cell lysates
 Lysates/proteins at 20 µg per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 87 kDa Observed band size: 75 kDa
 Exposure time: 4 seconds



All lanes: PLAP Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type A549 cell lysate
 Lane 2 : PLAP Rabbit mAb knockdown A549 cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human lung cancer tissue stained for PLAP using 60031 at 1/100 dilution in immunohistochemical analysis.



Immunocytochemistry/ Immunofluorescence PLAP antibody (60031)
 ICC/IF staining of PLAP in HeLa cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100.

Samples were incubated with 60031 at a working dilution of 1/100. The secondary antibody was Alexa FluorB 488 goat anti rabbit, used at a dilution of 1/500.

Nuclei

were counterstained with DAPI.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.