

ALKBH1 Rabbit mAb

Catalog No: #60064



Package Size: #60064-1 50ul #60064-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	ALKBH1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Human Mouse
Specificity	ALKBH1 Antibody detects endogenous levels of total ALKBH1
Conjugates	Unconjugated
Other Names	ABH; ABH1; alkB; ALKBH; ALKBH1; hABH;
Accession No.	Q13686
Calculated MW	44 kDa
SDS-PAGE MW	44 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200


Images



All lanes: ALKBH1 Rabbit mAb at 1/1k dilution
Lane 1 : JK whole cell lysates
Lane 2 : HepG2 whole cell lysates
Lane 3 : Mouse brain lysates
Lane 4 : Mouse liver lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 44 kDa
Observed band size: 44 kDa
Exposure time: 7 seconds



All lanes:ALKBH1 Rabbit mAb at 1/1k dilution
Lane 1 : Wild-type Hela cell lysate
Lane 2 :ALKBH1 Rabbit mAb knockdown Hela cell lysate
Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human colon tissue stained for ALKBH1 using 60064 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.