

MUTA Rabbit mAb

Catalog No: #60075



Package Size: #60075-1 50ul #60075-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	MUTA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Human Mouse Rat
Specificity	MUTA Antibody detects endogenous levels of total MUTA
Conjugates	Unconjugated
Other Names	MCM; Methylmalonyl CoA isomerase; Mut;
Accession No.	P22033
Calculated MW	83 kDa
SDS-PAGE MW	78 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200


Images




All lanes: MUTA Rabbit mAb at 1/1k dilution
Lane 1 : JK whole cell lysates
Lane 2 : HepG2 whole cell lysates
Lane 3 : MCF-7 whole cell lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 83 kDa
Observed band size: 78 kDa
Exposure time: 4 seconds



All lanes: MUTA Rabbit mAb at 1/1k dilution
 Lane 1 : Mouse heart lysates
 Lane 2 : Mouse liver lysates
 Lane 3 : Rat heart lysates
 Lane 4 : Rat liver lysates
 Lysates/proteins at 20 µg per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 83 kDa
 Observed band size: 78 kDa
 Exposure time: 6 seconds



All lanes: MUTA Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HeLa cell lysate
 Lane 2 : MUTA Rabbit mAb knockdown HeLa cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human cholangiocarcinoma tissue stained for MUTA using 60075 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.