

GRSF1 Rabbit mAb

Catalog No: #60107



Package Size: #60107-1 50ul #60107-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GRSF1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC;ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	GRSF1 Antibody detects endogenous levels of total GRSF1
Conjugates	Unconjugated
Other Names	G rich RNA sequence binding factor 1; G rich sequence binding factor 1; GRSF1;
Accession No.	Q12849
Calculated MW	53 kDa
SDS-PAGE MW	53 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500-1:2000


IHC: 1:50-1:200

ICC/IF: 1:50-1:200


Images




All lanes: GRSF1 Rabbit mAb at 1/1k dilution
Lane 1 : 293 whole cell lysates
Lane 2 : JK whole cell lysates
Lane 3 : Mouse brain lysates
Lane 4 : Rat brain lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 53 kDa
Observed band size: 53 kDa
Exposure time: 5 seconds



All lanes: GRSF1 Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HeLa cell lysate
 Lane 2 : GRSF1 Rabbit mAb knockdown HeLa cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human esophagus cancer tissue stained for GRSF1 using 60107 at 1/100 dilution in immunohistochemical analysis.



Immunocytochemistry/ Immunofluorescence GRSF1 antibody (60107)
 ICC/IF staining of GRSF1 in CEM cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100.

Samples were incubated with 60107 at a working dilution of 1/100. The secondary antibody was Alexa Fluor 488 goat anti rabbit, used at a dilution of 1/500.

Nuclei were counterstained with DAPI.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.