

Cip4 Rabbit mAb

Catalog No: #60117



Package Size: #60117-1 50ul #60117-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	Cip4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Human Mouse
Specificity	Cip4 Antibody detects endogenous levels of total Cip4
Conjugates	Unconjugated
Other Names	Cdc42 interacting protein 4; CIP4; DCIP4; hSTP; STOT; STP; Thyroid receptor interacting protein 10; trip10;
Accession No.	Q15642
Calculated MW	68 kDa
SDS-PAGE MW	75 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200


Images



All lanes: Cip4 Rabbit mAb at 1/1k dilution
Lane 1 : HepG2 whole cell lysates
Lane 2 : 3T3 whole cell lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 68 kDa
Observed band size: 75 kDa
Exposure time: 8 seconds



All lanes: Cip4 Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HAP1 cell lysate
 Lane 2 : Cip4 Rabbit mAb knockdown HAP1 cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human breast cancer tissue stained for Cip4 using 60117 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.