

Mras Rabbit mAb

Catalog No: #60172

Package Size: #60172-1 50ul #60172-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Mras Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Mras Antibody detects endogenous levels of total Mras
Conjugates	Unconjugated
Other Names	Mras; R ras3; Ras related protein MRas; Ras related protein RRas3; RRas3; Xras;
Accession No.	O14807
Calculated MW	24 kDa
SDS-PAGE MW	24 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500-1:2000

Images



All lanes: Mras Rabbit mAb at 1/1k dilution

Lane 1 : Hela whole cell lysates

Lane 2 : Mouse brain lysates

Lane 3 : Rat brain lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 24 kDa

Observed band size: 24 kDa

Exposure time: 9 seconds



All lanes:Mras Rabbit mAb at 1/1k dilution
Lane 1 : Wild-type
HeLa cell lysate
Lane 2 :Mras Rabbit mAb knockdown HeLa
cell lysate
Lysates/proteins at 20 µg per lane.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.