

PPM1G Rabbit mAb

Catalog No: #60173



Package Size: #60173-1 50ul #60173-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	PPM1G Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Human Mouse Rat
Specificity	PPM1G Antibody detects endogenous levels of total PPM1G
Conjugates	Unconjugated
Other Names	PP2C; PP2C gamma; PP2CG; PPM1G; PPP2CG; Protein phosphatase 1C; Protein phosphatase 1G;
Accession No.	O15355
Calculated MW	59 kDa
SDS-PAGE MW	73 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200


Images



All lanes: PPM1G Rabbit mAb at 1/1k dilution
Lane 1 : Hela whole cell lysates Lysates/proteins at 20 µg per lane.
SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 59 kDa Observed band size: 73 kDa
Exposure time: 5 seconds



All lanes: PPM1G Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HeLa cell lysate
 Lane 2 : PPM1G Rabbit mAb knockdown HeLa cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human colon carcinoma tissue stained for PPM1G using 60173 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.