

MDH2 Rabbit mAb

Catalog No: #60180



Package Size: #60180-1 50ul #60180-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	MDH2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC
Species Reactivity	Human Mouse Rat
Specificity	MDH2 Antibody detects endogenous levels of total MDH2
Conjugates	Unconjugated
Other Names	M MDH; Malate dehydrogenase;
Accession No.	P40926
Calculated MW	36 kDa
SDS-PAGE MW	36 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200


Images



All lanes: MDH2 Rabbit mAb at 1/1k dilution
Lane 1 : K562 whole cell lysates Lysates/proteins at 20 µg per lane.
SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 36 kDa Observed band size: 36 kDa
Exposure time: 5 seconds



All lanes:MDH2 Rabbit mAb at 1/1k dilutionLane 1 : Wild-type HAP1 cell lysateLane 2 :MDH2 Rabbit mAb knockdown HAP1 cell lysateLysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human breast cancer tissue stained for MDH2 using 60180 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.