

## Ubiquitin K48 Rabbit mAb

Catalog No: #60303

Package Size: #60303-1 50ul #60303-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Ubiquitin K48 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG,Kappa
Purification	Protein A
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Endogenous
Conjugates	Unconjugated
Target Name	Ubiquitin K48
Modification	Unmodified
Other Names	UBB;Polyubiquitin-B;UBC;Polyubiquitin-C;RPS27A;UBA80;UBCEP1;Ubiquitin-40S ribosomal protein S27a;Ubiquitin carboxyl extension protein 80;UBA52;UBCEP2;Ubiquitin-60S ribosomal protein L40;CEP52;Ubiquitin A-52 residue ribosomal pr;Ubiquitin K48;
Accession No.	P0CG47;P0CG48;P62979;P62987
Calculated MW	26kDa
SDS-PAGE MW	13-250kDa
Formulation	PBS; 50% glycerol; 0.05% Proclin 300; 0.05%BSA
Storage	-15°C to -25°C/1 year(Do not lower than -25°C)

## Application Details

IHC 1:200-1:1000;WB 1:500-1:2000;IF 1:200-1:1000;ELISA 1:5000-1:20000;

## Images



Western blot analysis of lysates from HEK293 cells, using Ub Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.