

ApoE Rabbit mAb

Catalog No: #60917



Package Size: #60917-1 50ul #60917-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ApoE Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG, Kappa
Purification	Recombinant Antibody expressed in animal component-free (ACF) media, purified via Protein A affinity chromatography.
Applications	WB;IHC;IF;IP;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Endogenous
Immunogen Description	The specific immunogen used to produce this antibody is proprietary information.
Conjugates	Unconjugated
Target Name	APOE
Other Names	APOE ; Apolipoprotein E ; Apo-E
Accession No.	P02649
Calculated MW	36kDa
SDS-PAGE MW	36kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year


Application Details

WB 1:2000-1:10000; IHC 1:500-1:2000; IF 1:200-1:1000; IP 1:50-1:200; ELISA 1:5000-1:20000

Images



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-ApoE antibody.



Human liver was stained with anti-ApoE Rabbit antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.