

NGAL Rabbit mAb

Catalog No: #60918



Package Size: #60918-1 50ul #60918-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NGAL Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG, Kappa
Purification	Recombinant Antibody expressed in animal component-free (ACF) media, purified via Protein A affinity chromatography.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Endogenous
Conjugates	Unconjugated
Target Name	LCN2 HNL NGAL
Other Names	Neutrophil gelatinase-associated lipocalin ; NGAL ; 25 kDa alpha-2-microglobulin-related subunit of MMP-9 ; Lipocalin-2 ; Oncogene 24p3 ; Siderocalin LCN2 ; p25 ;
Accession No.	P80188
Calculated MW	23kDa
SDS-PAGE MW	23kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year


Application Details

WB 1:2000-1:10000; IHC 1:1000-1:4000; IF 1:200-1:1000; ELISA 1:5000-1:20000

Images



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-NGAL antibody.



Human spleen was stained with anti-NGAL Rabbit antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.