

C3 Rabbit mAb

Catalog No: #60919

Package Size: #60919-1 50ul #60919-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG, Kappa
Purification	Recombinant Antibody expressed in animal component-free (ACF) media, purified via Protein A affinity chromatography.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Endogenous
Conjugates	Unconjugated
Target Name	C3
Other Names	CPAMD1
Accession No.	P01024
Calculated MW	187kDa
SDS-PAGE MW	187kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year


Application Details

WB 1:1000-1:5000; IHC 1:200-1:1000; IF 1:200-1:1000; ELISA 1:5000-1:20000


Images




Various whole cell lysates were separated by 4-20% SDS-PAGE, and the primary antibody was used at 4°C over night with a 1:2500 dilution .



Human brain was stained with anti-C3 Rabbit antibody



Mouse brain was stained with anti-C3 Rabbit antibody



Rat brain was stained with anti-C3 Rabbit antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.