

Anti-Human CD366/HAVCR2/TIM-3 Recombinant Nanobody

Catalog No: #60922

Package Size: #60922 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-Human CD366/HAVCR2/TIM-3 Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Purification	Purified by Nickel column
Applications	ELISA, FCM
Species Reactivity	Human
Target Name	T-cell immunoglobulin mucin receptor 3, HAVcr-2, TIMD-3, T-cell membrane protein 3, HAVCR2, Hepatitis A virus cellular receptor 2, CD366, TIMD3, TIM-3, TIM3, T-cell immunoglobulin and mucin domain-containing protein 3
Accession No.	Q8TDQ0
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.ELISA 1:2000-1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.