

TFE3 Mouse mAb

Catalog No: #64443

Package Size: #64443-1 50ul #64443-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TFE3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1, Kappa
Purification	Recombinant Antibody expressed in animal component-free (ACF) media, purified via Protein G affinity chromatography.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Bovine
Specificity	This antibody detects endogenous levels of TFE3 protein.
Immunogen Description	Synthesized peptide derived from human TFE3
Conjugates	Unconjugated
Target Name	TFE3 BHLHE33
Other Names	TFE3 ; BHLHE33 ; Transcription factor E3 ; Class E basic helix-loop-helix protein 33 ; bHLHe33
Accession No.	P19532
Calculated MW	62kDa
SDS-PAGE MW	62kDa
Concentration	1 mg/ml
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year


Application Details

WB 1:500-2000;IHC 1:200-1000;IF 1:100-500;ELISA 1:1000-5000

Images



Various whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-TFE3 antibody.



Human testis tissue was stained with Anti-TFE3 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.