

alpha Synuclein Antibody Biotin Conjugated

Catalog No: #C00005B

Package Size: #C00005B 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	alpha Synuclein Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 45-95 140 derived from human alpha Synuclein
Conjugates	Biotin
Target Name	alpha Synuclein
Other Names	PD1; NACP; PARK1; PARK4; Alpha-synuclein; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; SNCA
Cell Localization	Intracellular
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Immunohistochemistry1:100-500

Background

May be involved in the regulation of dopamine release and transport. Induces fibrillization of microtubule-associated protein tau. Reduces neuronal responsiveness to various apoptotic stimuli, leading to a decreased caspase-3 activation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.