

Neurokinin A Antibody FITC Conjugated

Catalog No: #C00051F

Package Size: #C00051F 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurokinin A Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 98-107 129 derived from human NKA
Conjugates	FITC
Target Name	Neurokinin A
Other Names	NK2; NPK; NKNA; TAC2; Hs.2563; Protachykinin-1; PPT; TAC1; NKA
Accession No.	Swiss-Prot#P20366NCBI Gene ID6863
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Tachykinins are active peptides which excite neurons, evoke behavioral responses, are potent vasodilators and secretagogues, and contract (directly or indirectly) many smooth muscles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.